

Current Status And Results Of The Terumo Relay Branched Endograft Device For Aortic Arch Lesion Repairs: Unique Features And When Will It Be available In The US

TIGER Post Market Registry

V. Rimbau, MD, PhD, RAMC
On behalf of the TIGER Collaborators
Prof. and Chief of Vascular Surgery Department and Teaching Vascular Unit
Cardiovascular Institute, Hospital Clinic of Barcelona
University of Barcelona





Disclosures

I have the following potential conflicts of interest to report:

- Consulting: Terumo Aortic, Medtronic, Cook Medical, iVascular, W.L. Gore, VB Devices, Lifetech, Möller Medical, Aortyx, Balt Medical.
- Stockholder of a healthcare company: Aortyx, VB Devices
- PI of TIGER Postmarket Registry






TIGER Global Site Overview




As of April 4th, 2024
Total Subject Enrollment is 1190

52 Sites active across 10 Countries

Site Enrollment by County

- o Spain - 4 sites
- o Portugal - 2 sites
- o Netherlands - 7 sites
- o Germany - 7 sites
- o Belgium - 1 sites
- o France - 2 sites
- o Italy - 7 sites
- o Switzerland- 1
- o UK - 3 sites
- o USA - 18 sites








TIGER Overview

- **Study Size**
 - Minimum of **1000 patients** at up to **80 sites** globally
 - No upper limit set
- **Study Duration**
 - **Enrolment opened Nov 2019**, no closure date
 - Data will be collected at the following time points:
 - o Baseline (Procedure)
 - o Early follow-up (Discharge/30 days to six months)
 - o One-year post index procedure
 - o Post index procedure visit: 1st long-term follow-up annually through subject lifetime or until lost to follow up (LTFU)
 - For Thoracic cohort, a **minimum of 5 Years** annual follow-up is required

Key Study Contacts


- Prof Vincent Rimbau – Hospital Clinic i provincial de Barcelona (Coordinating Investigator – Global)
- Prof Michele Reijnen -Rijnstate, Arnhem, Netherlands (Regional PI – Europe)
- Mr David Murray – Manchester University Hospital, United Kingdom (CI - UK)








Module Overview

One overarching protocol utilising a modular



- All patients included
- Common Endpoints

• FDA Post Approval Studies for Aneurysm & Dissection Indications for RelayPro embedded into Registry








Relay® Branch Thoracic Stent-Graft System (Outside US)

- Builds upon the **core features** of the standard Relay®Pro platform
- Incorporating the same advance technology and design principles
- Added customization*
- Benefits of this platform include

- Graft**
 - Reinforced endograft material
 - Reinforced endograft
 - Low cover technology (where necessary)
- Delivery System**
 - Fast track technology
 - Pre curved distal inner catheter
 - Repositionable
 - Experimental proximal clipping
- Branches**
 - With reinforced design
 - Proximal clipping and distal open core delivery system
 - Customizable design
 - US FDA delivery system

- Delivery system pre-curved inner catheter**
 - Experimental proximal clipping
- Dual sheath technology**
 - Experimental proximal clipping
- Precise proximal landing and progressive apposition**
 - Experimental proximal clipping

Procedural (Site Reported)

	% (n/N) or Mean ± SD (n)	Range
Surgical Acuity		
• Elective	100% (19/19)	
• Emergent	0	
Anesthesia Type		
• General Anesthesia	100% (19/19)	
• Local Anesthesia	0	
• Other	0	
Mean Hospital Stay (days)	12.5 ± 11.7 (19)	2 – 52
Mean ICU Stay (hours)	88.3 ± 97.9 (19)	72 – 336
Total Surgery Time (min)	292.8 ± 116.3 (19)	268 – 540
Total Endovascular Time (min)	206.5 ± 73.8 (17)	180 – 360
Estimated Blood Loss (mL)	772.2 ± 1028.7 (9)	500 – 3500
Transfusion Required	36.8% (7/19)	
Vascular Access		
• Cut Down	36.8% (7/19)	
• Percutaneous	47.4% (9/19)	
• Conduit	15.8% (3/19)	
Discharge Destination		
• Temporary Ambulatory Care Center	0	
• Return to pre-op living location	89.5% (17/19)	
• Other	10.5% (2/19)	

Technical Success (Site Reported)

Components of Technical Success	% (n/N)
• Device introduction into the entry site acceptable	100% (18/18)
• Advancement of device to the lesion site acceptable	100% (18/18)
• Device deployment from the delivery system acceptable	100.0% (18/18)
• Accuracy of the device system deployment acceptable	100.0% (18/18)
• Stent-graft patent	100.0% (18/18)
• Device integrity maintained (no tears, fractures, etc.)	100% (18/18)
• Stent-graft deployed without kinking or twisting	100% (18/18)
• Absence of Type Ia, Type Ib, Type IIIa and Type IIIb endoleak	94.7% (18/19)
• Absence of unplanned coverage of aortic branch vessels	100% (19/19)

100% Acceptable Device Introduction, Deployment, Accuracy, Patency, and Integrity

Mortality (Site Reported)

	Procedure (POD 0 – POD 29)	Early Follow-up	12 Months	2 Years	3 Years	Total (subjects)
Visits Performed	19	8	4	3	1	19
All-Cause Mortality	2	1	0	0	0	3
Aortic-Related Mortality	2	0	0	0	0	2

No Aortic-related mortality was reported as device related

- Aortic Related Death (n=2)
- Definition: Aortic-related death is defined as any death occurring within 30 days of implant, due to rupture or following any procedure intended to treat the target lesion
 - Subject 004-016 (POD 13), acute liver failure
 - Subject 004-018 (POD 9), heart failure

Stroke (Site Reported)

	After Index Procedure (POD 1 – POD 29)	Early Follow-up	12 Months	2 Years	3 Years
Visits Performed	19	8	4	3	1
Stroke	0	0	0	2	0

Stroke (n=2)


- Subject 015-008, 2-year visit (01/27/2022), Site reported slight disability. Able to look after own affairs without assistance, but unable to carry out all previous activities
- Subject 024-028, 2-year visit (01/29/2024), Site reported a moderately severe disability. Unable to attend to own bodily needs without assistance & unable to walk unassisted

Secondary Interventions (Site Reported)

	After Index Procedure (POD 1 – POD 29)	Early Follow-up	12 Months	2 Years	Total
Visits Performed	19	8	4	3	19
Subjects w/ any Intervention	2	0	2	0	4
Total number of Interventions	2	0	2	0	4




Reasons for Secondary Interventions (Site Reported)

- 4 subjects w/ 4 Secondary Interventions at Follow-up**
- **Major bleeding (POD 1)** - Wound revision to achieve hemostasis of the arterial access
- **Type II Endoleak and Stenosis (POD 3)** - Reintervention performed for Type II endoleak
- **Proximal kinking and suspected Type I Endoleak (POD 259)** - embolization of the false lumen of the dissection, and dilatation of the slight collapse of the proximal edge of the endoprosthesis were performed
- **Extension of previous TEVAR (POD 275)** - Extension of previous TEVAR due to distal aneurysm with an increased diameter of 47 to 55 mm



More data...

- No ruptures or conversions to open surgery
- No device occlusions
- No clinically significant migration >10mm
- No wire frame fractures








Summary

A total of **19 subjects** are enrolled in TIGER Relay@Branch

- ❖ Enrollment and follow-up remain ongoing across all TIGER Modules
- ❖ Most of the components of technical success (e.g., deployment, accuracy, delivery, etc.) were 100%
- ❖ No aortic-related mortalities related to the device or lesion
- ❖ No ruptures or conversions to open surgery
- ❖ No stent fractures or migration >10mm

Preliminary results continue to support the safety and effectiveness of the Relay@Branch Thoracic Stent-Graft System



Current Status And Results Of The Terumo Relay Branched Endograft Device For Aortic Arch Lesion Repairs: Unique Features And When Will It Be available In The US

TIGER Post Market Registry

V. Riambau, MD, PhD, RAMC
On behalf of the TIGER Collaborators

Prof. and Chief of Vascular Surgery Department and Teaching Vascular Unit
Cardiovascular Institute, Hospital Clinic of Barcelona
University of Barcelona



