




One Year Results Of The SPORTs Trial: DCBs Vs. DES (Eluvia) VS. Bare Stents In Long SFA Lesions (TASC C AND D) - What Is Best?

Prof. Dr. med. Gunnar Tepe,
Institut für Diagnostische und Interventionelle Radiologie
RoMed Klinikum Rosenheim



Disclosure

Study support by B.Braun, Biotronic, BSC, CVT, Philips, Medtronic, Shockwave




Within the prior 24 months, I have had a relevant financial relationship(s) with an ineligible company(ies) listed below.



Nature of Financial Relationship	Ineligible Company
Grant/Research Support	BSC, B Braun, Biotronic, CVT, Medtronic, Medalliance, Philips
Consultant Fees/Honoraria	BSC, Medtronic
Individual Stock(s)/Stock Options	
Royalties/Patent Beneficiary	
Executive Role/Ownership Interest	
Other Financial Benefit	




Limitation on the current data TASC C + D for endovascular intervention



- Most prospective randomized studies
 - IDE trials, TASC A and B
- Truly long lesions
 - Often retrospective or registries
 - No control group and dubious indirect comparisons
 - No standardized f/u, no core lab

No definitive data on primary strategy decision making





Current treatment possibilities for TASC C + D

- Endovascular
 - POBA
 - DCB (+ high rate bail out stent)
 - Stent (which one?)
 - DES
 - Covered Stentgraft
- Bypass Surgery






SPORTS



- DESIGN: Prospective, randomized, open-label, 3-arm, multi-center
- OBJECTIVE: To compare angiographic and clinical outcome of TASC C and D lesions in the SFA after treatment with PCB, PES or BNS
- PRIMARY ENDPOINT: Angiographic diameter stenosis at 12 month
- HYPOTHESES: 1) Superiority of PES over BNS
2) Non-inferiority of PCB versus BNS
- PRINCIPAL INVESTIGATOR
Gunnar Tepe, MD
RoMed Klinikum, Germany

PCB: paclitaxel coated balloon (Siemens Please
OTW);
PES: paclitaxel eluting VEITH AVE stent
system);
BNS: Biotronic Biotricor stent system



SPORTS

MAIN INCLUSION CRITERIA


- Rutherford classes 2 – 4
- Peripheral lesions in the SFA and / or A.pop. seg I
- Minimum diameter stenosis $\geq 70\%$
- Treatment length at least 15 cm (lesion length at least 13 cm)

MAIN EXCLUSION CRITERIA

- ISR lesions
- A.pop. Seg II or below or within a bypass graft
- Acute thrombosis of the study lesion
- Strongly calcified lesions with circumferential presence of calcifications and a lesion calcification length >4 cm
- RVD <4 mm and >6 mm
- Lesions below the knee requiring treatment
- No patent distal run-off vesse!

SPORTS



Prof. S. Müller-Hulsbeck
DIAKO Krankenhaus Flensburg

Prof. J.P. Goltz
SANA Klinikum, Lübeck

Prof. J.P. Goltz/Dr. Stroth
Universitätsklinikum Schleswig-Holstein, Lübeck

Dr. M. Lichtenberg
Klinikum Hochsauerland, Arnsberg

Prof. T. Zeller
Universitäts-Herzzentrum Bad Krozingen




Dr. H. Schröder
Ihre-Radiologen, Jüdisches Krankenhaus, Berlin

Prof. K. Brechtel
Ihre-Radiologen, Franziskus Krankenhaus, Berlin

Prof. T. Abrecht
Vivantes Klinikum Neukölln, Berlin

Prof. G. Tepe
RoMed Klinikum Rosenheim



Prof. H. Brodmann
Medizinische Universität Graz

Patient baseline characteristics



	PCB n = 74	PES n = 74	BNS n = 76
Male (%)	66	60	72
Age (years, mean±SD)	70 ± 9	68 ± 8	67 ± 9
Weight (kg, mean±SD)	78 ± 19	79 ± 16	79 ± 16
Diabetes (%)	30	23	26
Hypertension (%)	88	80	78
Hyperlipidemia (%)	66	70	71
Obesity (%)	18	12	21
Renal disease (%)	12	8	3
Current smoker (%)	60	55	58
ABI study leg (mean±SD)*	0.6 ± 0.2 n = 60	0.6 ± 0.2 n = 64	0.6 ± 0.2 n = 61

RUTHERFORD CLASSIFICATION BASELINE



Lesion baseline characteristics

	PCB n = 74	PES n = 74	BNS n = 76	P value
Lesion length, mm, mean ± SD	221 ± 87 n=73	235 ± 78 n=73	227 ± 78 n=75	0.57
Occlusion, n (%)	52 (70)	63 (85)	56 (74)	0.08
Occlusion length, mm, mean ± SD	175 ± 91 n=51	179 ± 89 n=62	151 ± 81 n=54	0.18
RVD, mm, mean ± SD	5.0 ± 0.6	5.3 ± 0.7	5.2 ± 0.7 n=75	0.01
MLD in lesion, mm, mean ± SD	0.4 ± 0.7	0.2 ± 0.6	0.3 ± 0.7	0.18
Diameter stenosis in lesion, %, mean ± SD	92.6 ± 13.2	96.8 ± 9.7	94.2 ± 11.7	0.10

Lesion baseline characteristics

	PCB n = 74	PES n = 74	BNS n = 76	P value
Calcification (PACCS), n (%) - Corelab				
0	3 (4.1)	2 (2.7)	5 (6.6)	0.36
1	15 (20.3)	28 (37.8)	18 (23.7)	
2	1 (1.4)	0	0	
3	38 (51.4)	29 (39.2)	31 (40.8)	
4	15 (20.3)	14 (18.9)	20 (26.3)	
NA	2 (2.7)	1 (1.4)	2 (2.6)	

Procedural data

	PCB n = 74	PES n = 74	BNS n = 76	P value
Pre-dilation (%) ¹	91	95	90	0.50
Ball out stenting (%) ¹	58	NA	NA	NA
Treated length, mm (mean±SD) ²	269 ± 82 n=73	273 ± 80 n=73	257 ± 68 n=74	0.42
Post-procedural MLD in lesion, mm (mean±SD) ²	3.5 ± 0.7 n=73	4.2 ± 0.7 n=74	4.1 ± 0.7 n=74	<0.0001
Post-procedural diameter stenosis, % (mean±SD) ²	30 ± 12 n=73	20 ± 15 n=74	20 ± 14 n=74	<0.0001

¹ CRF data; ² Corelab analysis; NA = not applicable

Efficacy and safety endpoints

12 month Corelab data, ITT



	PCB n = 48	PES n = 56	BNS n = 52	P value
Diameter stenosis in lesion, %, mean ± SD	53.7 ± 33.9	25.3 ± 32.4	60.0 ± 29.8	
LLL in lesion, mm, mean ± SD	1.1 ± 1.8	0.4 ± 1.7	2.0 ± 1.5 n=50	<0.0001
MLD in lesion, mm, mean ± SD	2.3 ± 1.7	3.9 ± 1.7	2.1 ± 1.6	<0.0001

Hypothesis testing for primary endpoint (in-lesion diameter stenosis)

Hypothesis 1: Superiority DES vs. BNS ✔ $p=0.0001$ (-34.7 [CI 97.5%: -47.7]; -20.2)¹*

Hypothesis 2: Non-inferiority PCB vs. BNS ✔ $p=0.0010$ (-6.3 [CI 97.5%: -22.3; 7.6])¹*

*Intention-to-treat and p-values based on per-protocol analysis, no covariate



- ### Summary and Conclusion
- PES showed significantly better angiographic results 12 months after treatment compared to treatment with BNS
 - Non-inferiority of PCB treatment compared to BNS shown for angiographic stenosis degree
 - PCB treatment in long, complex lesions required high bailout stenting rate, mainly due to dissections
 - TLR rates favor PES over PCB and BNS treatment
 - Long-term data necessary
 - Subgroup analyses planned