

The Role Of Pressure Wires And IVUS In The Assessment Of Renal Artery Stenosis In The Treatment Of Hypertension And Other Conditions

Lindsay Machan, MD
University of British Columbia
Vancouver, BC, Canada



Disclosures

- Medical Advisory Board
 - Auxetics
 - Boston Scientific
 - Medtronic
- Stock Shareholder
 - Analytics 4 Life
 - Auxetics
 - Calgary Scientific
 - Cohesic
 - Harmonic Medical
 - Ikomed Medical
 - Innovere Medical
 - Zymeworks

Prediction of Hypertension Improvement After Stenting of Renal Artery Stenosis
Comparative Accuracy of Translesional Pressure Gradients, Intravascular Ultrasound, and Angiography


Massoud A. Lecar, MD,* Jai Varma, MD,* Adam Shapiro, MD,* Ibrahim Fahash, MD,* Seyed T. Raza, MD,† Ziad Elghoul, MD,* Anthony C. Leonard, PhD,‡ Karthikeyan Meganathan, MS,§ Sohail Ikram, MD*

- 62 patients
- Compared:
 - Pressure guidewire
 - FFR
 - Mean gradient
 - Hyperemic systolic gradient (after 30 mg IA papaverine)
 - IVUS
 - Quantitative angiography

JACC. 53: 2009;2363-71

Prediction of Hypertension Improvement After Stenting of Renal Artery Stenosis
Comparative Accuracy of Translesional Pressure Gradients, Intravascular Ultrasound, and Angiography


Massoud A. Lekar, MD,* Jai Varma, MD,* Adam Shapiro, MD,* Ibrahim Fahash, MD,* Seyed T. Raza, MD,† Ziad Elghoul, MD,* Anthony C. Leonard, PhD,‡ Karthikeyan Meganathan, MS,§ Sohail Ikram, MD*




- HSG ≥ 21 mm Hg only independent predictor HBP improvement
 - (odds ratio: 1.39; 95% confidence interval: 1.05 to 1.65; p 0.013)

JACC. 53: 2009;2363-71


39-year-old male with medically refractory hypertension and right renal artery stenosis



39-year-old male with medically refractory hypertension and right renal artery stenosis




39-year-old male with medically refractory hypertension and right RAS



Angioplasty - 6 x 20 mm balloon No change angio or pressures

39-year-old male with medically refractory hypertension and right RAS



Angioplasty - 6 x 20 mm balloon No change angio or pressures 7 x 18 mm BeGraft covered stent

39-year-old male with medically refractory hypertension and right RAS




Post BeGraft

39-year-old male with medically refractory hypertension and right RAS



Post BeGraft




IVUS - stent 5.8 mm


39-year-old male with medically refractory hypertension and right RAS



Post BeGraft




IVUS - stent 5.8 mm



Post dilated to 7.0 mm
Distal to stent: 143/74, mean 99
Aorta: 135/79, mean 99


84-year-old male solitary right kidney

- Progression CKD
 - Cr 495, GFR 9
- ER visits
 - SOB due to fluid overload
 - Worsening hypertension



Ultrasound impossible and MR poor due to dyspnea


84-year-old male solitary right kidney



Ultrasound impossible and MR poor due to dyspnea

Pressure wire no gradient


48 yr old woman
HBP poorly controlled on 3 medications. CVA x 2



FMD right renal artery

LRA probably normal


48 yr old woman
HBP poorly controlled on 3 medications. CVA x 2



Angio -> FMD

Hilum: 161/93, mean 120 mm Hg
Aorta: 129/80, mean 101

48 yr old woman
HBP poorly controlled on 3 medications. CVA x 2




Angio -> FMD

Hilum: 161/93, mean 120 mm Hg

Balloon angioplasty

48 yr old woman
HBP poorly controlled on 3 medications. CVA x 2




Angio -> FMD

Hilum: 161/93, mean 120 mm Hg
Aorta: 129/80, mean 101

Balloon angioplasty

Post dilation angio
Persistent undulation - no gradient

48 yr old woman
HBP poorly controlled on 3 medications. CVA x 2



Angio -> subtle FMD

18 mm systolic gradient

Balloon angioplasty

31 yr old man. Type A dissection. Hypertension on five medications.
Query dynamic dissection flap across renal artery orifice

Chronic occlusive changes RRA
Pressure wire – no gradient
IVUS - no **dynamic** dissection flap

Aortic Pullback Pressures

- 31 mm Hg systolic gradual gradient through narrowed true lumen
- Patient elected medical management

Summary - The Role Of Pressure Wires And IVUS In The Assessment Of Renal Artery Stenosis In The Treatment Of Hypertension And Other Conditions

- Pressure wires
 - Guide intra-procedure management atherosclerotic RAS
 - ? More accurate than IVUS threshold for intervention
 - FMD – direct intervention in subtle cases
 - Good endpoint when angio findings unchanged
- IVUS
 - Accurate sizing endovascular devices
 - Dynamic assessment, e.g. dynamic flaps dissection