

Update on F/BEVAR for TAAAs degenerative and PDTAAAs: The Value of Inner Branched Devices and the Be-Back Catheter to Cross the Dissection Flap and for other Purposes

E.L.G. Verhoeven, MD, PhD, and A. Katsargyris, MD
Department of Vascular and Endovascular Surgery, Paracelsus Medical University, Nuremberg, Germany


Disclosures

- William Cook Europe/Cook Inc.
 - Consultant & Research grants
- Bentley
 - Consultant

Degenerative TAAA


- Basic Rule:
 - Find/Create long and durable proximal landing zone
- Anatomical Difficulties
 - Quality of target vessels
 - Angulations
- Current advantages:
 - Upper and Femoral approach for Branches
 - Multitude of BE and SE Stents for bridging
 - Inner Branches (and multiple graft manufacturers)

F/B-Graft with both Fens and Branches Customize the design to the patients' Anatomy






Inner Branch (as a third option)

- When
 - Limited space between graft and aortic wall
 - Target vessels with sharp take-off
- Combination Fenestrations/Inner Branches
 - N=83
- Inner branches only
 - N=6



Additional Risk due to sharper Angle? (Especially in Renal Arteries?)





Journal Pre-proof

Midterm single center results with the use of custom-made endografts with inner branches, a call for attention.



Jose I. Torrealba, MD, Giuseppe Panuccio, MD, PhD, Petroula Nana, PhD, Antonino Giordano, MD, Daour Yousef Al Sarhan, MD, Tilo Köbel, MD, PhD

Key Findings: Custom-made inner branch devices were used in 69 patients with a high bridging stent occlusion rate, especially in renal arteries (17% at 12 months). Misalignment of the inner branch opening, and preoperative acute aortic angles were identified as related factors leading to occlusions.

- ### PDTAAA
- #### Specific Anatomical Features with regard to F/BEVAR (More forgiving than we thought!)
- Access
 - Small True Lumen
 - Target Vessels originating from False lumen
 - Soft and Wide Target Vessels
 - Distal Landing Zone in dissected Artery

- ### PDTAAA
- #### Specific Anatomical Features with regard to F/BEVAR (More forgiving than we thought!)
- Access
 - Small True Lumen
 - Target Vessels originating from False lumen
 - Soft and Wide Target Vessels
 - Distal Landing Zone in dissected Artery

- ### Target Vessels from False Lumen Perforating the dissection flap?
- (Steerable) sheaths
 - Back of a wire
 - TIPPS needle
 - BeBack Re-entry Catheter
-



Be-Back Re-entry Catheter



- Steerable Needle/Flexible Catheter
- Re-entry:
 - Perforation of dissection flap in PDTAAA
 - In-situ fenestration LSA (in acute cases)
- Tricks:
 - Support with Steerable/guiding sheath
 - In situ: balloon inside graft





BeBack crossing catheter
The unique design with tapered re-entry capabilities

3.0x37 model with 0.014" guidewire lumen

BeBack

2 needle lengths and applications (2.9 F)





- Conclusions
- F/B grafts for TAAA and PDTAAA
 - Customize graft to the anatomy/Off-the-shelf where possible/needed
 - ➔ Do not be dogmatic
 - Inner branches valuable third option
 - ➔ In selected Vessels only!
 - Be-Back steerable needle
 - ➔ Very useful to perforate dissection flap & for acute ISF